


TOTAL INTERNAL REFLECTION
6A44.20



SET UP ON A LARGE OPTICAL BENCH ON A CART.

USE He-Ne LASER AND GREEN ACRYLIC SHEET ON NORTH WALL.

LASER MUST BE VERY CLOSE TO ACRYLIC SHEET.

MAKE SURE THE ACRYLIC SHEET IS STABLE BEFORE MOVING.

END SUPPORTS ARE SQUARE U-SHAPED HOLDERS, MIDDLE SUPPORT IS A THREE-FINGERED SPRING CLAMP.